





MILLENNIUM SERIES - E022G1L

Progressive Cavity Pump



Millennium E022G1L



Features

- Low Shear
- Solids Handling up to 50%
- Hardened Joints
- Packing or Mechanical Seal
- 3-Stage Options
- Steel or Stainless Steel

Materials

Rotor: Chrome-Plated Tool Steel or Stainless Steel

Stator: Buna Nitrile, Natural Rubber, Viton™, EPDM

Motor: 1-1/2 - 10 hp

Inlet: 6"

Outlet: 6" - ANSI Flange

Bases: Steel and Stainless Steel

Pump Description

The Liberty Millennium Series E022G1L Heavy Duty Progressive Cavity Pump is ideal for high viscosity liquids and heavy waste solids and other abrasive pumping applications such as blood and waste in poultry processing facilities. The E022G1L's flow rate of 22 gallons/100 revolutions @ 55 to 99 GPM coupled to its high temperature resistance of up to 350° degrees Fahrenheit makes it an ideal solution for your industrial heavy-duty pumping needs.

The E022G1L is a standard ANSI flanged pump design manufactured with high strength cast iron or 316 Stainless Steel pump casings and designed in 1, 2 and up to 4 stages for 85, 170, 350 psi discharge pressures making it an ideal low-cost, maintenance free drop-in replacement pump. Standard suction and discharge flange size for the E022G1L is 6".

The Millennium Series low shear and gentle pumping action enables it to pump delicate fluids and not emulsify the liquid making it an ideal polymer pump or oil/water separation pump. Solids handling capability of up to up to 50% by weight allows even the toughest application to be pumped. The Millennium's modular design comes with a positively sealed, hardened gear joint assembly, hardened tool steel rotors or stainless steel rotors both with hard chrome plating for long-life pumping and a wide selection of stators in Buna Nitrile, Viton, EPDM and Natural Rubber.

All other wetted parts are either carbon steel or 316 stainless steel. The standard seal design is a set of packing with a Lantern Ring and flush connectors. Mechanical seal options are readily available.

Contact a Liberty Process Equipment Representative today for more information about the Millennium E022G1L Series of Progressive Cavity Pumps for blood and waste in poultry processing plants.

Curves

Click Links Below

1-Stage
1E022G1L


2-Stage
2E022G1L

Drawings

Click Links Below

1-Stage
1E022G1L

2-Stage
2E022G1L



P: 847.640.7867 F: 847.640.7855
2525 Clearbrook Drive Arlington Heights, IL 60005



www.libertyprocess.com



MILLENNIUM SERIES - 1E022GIL

Progressive Cavity Pump

Performance Curves



1 Stage Pump

70 Durometer (Performance based on water)

[BACK TO TOP](#)

P: 847.640.7867 F: 847.640.7855
2525 Clearbrook Drive Arlington Heights, IL 60005



www.libertyprocess.com



MILLENNIUM SERIES - 2E022GIL

Progressive Cavity Pump

Performance Curves



2 Stage Pump

70 Durometer (Performance based on water)

[BACK TO TOP](#)

P: 847.640.7867 F: 847.640.7855
2525 Clearbrook Drive Arlington Heights, IL 60005


www.libertyprocess.com



MILLENNIUM SERIES - 1E022GIL

Progressive Cavity Pump

Dimensional Drawings



(Pump Rotation is Counter Clockwise) ALL DIMN: Inches

Net weight

≈ lbs.

440

BACK TO TOP

1 Stage Pump

P: 847.640.7867 F: 847.640.7855
2525 Clearbrook Drive Arlington Heights, IL 60005


www.libertyprocess.com



MILLENNIUM SERIES - 2E022GIL

Progressive Cavity Pump

Dimensional Drawings



BACK TO TOP

(Pump Rotation is Counter Clockwise) ALL DIMN: Inches

Net weight

≈ lbs.

520

2 Stage Pump

P: 847.640.7867 F: 847.640.7855
2525 Clearbrook Drive Arlington Heights, IL 60005

www.libertyprocess.com